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PREPARATION

1. With the old or new prosthesis, a functional impression is made by the dentist.
2. Produce plaster model and embed directly into the cuvette.
3. Remove excess impression material.
4. Insulate teeth with Optosil (or similar); Unconditionally release occlusal surfaces.
5. Insulate and counter plaster and denture.
6. Brew the cuvette in a water bath and open.
7. Remove impression material.
8. Mill out the functional edge and inner surface of the prosthesis approx. 1mm.
9. Leave the edges of the prosthesis. Always frame and bond Prothesoft® with a hard plastic layer.
10. Isolate denture and plaster and allow to dry completely.
11. Wet or degrease the inner surfaces of the prosthesis with a monomer that is available in the laboratory.

PREPARATION OF THE GEL

1. Add the contents of bottle B (peroxide powder) to bottle A (basic substance monomer).
2. Close bottle A and shake well to mix both substances.
3. Keep lid of bottle A closed; Cover the lid with your hand,
4. Heat the bottle in a water bath (90 ° C) and shake well at intervals of approx. 7-10 seconds. Constantly check the consistency so as not to miss the time when the gel becomes viscous.
5. Within 1 ½ - 2 minutes, a thick gel is formed in the bottle (honey resembles).
6. Open the bottle, avoiding inflowing water! Pour gel on plaster model or denture and spread evenly with a spatula.
7. Close the cuvette and place under a press as usual (about 3 minutes). Tighten the clamp tightly.
8. Remove excess mass with dry paper.
9. Polymerize the cuvette in a water bath (100 ° C) for approx. 2 hours.
10. Then remove the cuvette from the water bath and allow to cool. Under no circumstances should this happen in cold water!
11. Open the cuvette, remove the denture and work out and polish as usual.

REMARKS

1. If the relining is still too soft after opening the cuvette, it must be dried in the air for about 30 minutes. Then it will be edited.
2. The polymerization time in the water bath must not be shorter than 90 minutes to get rid of the residual monomers.
3. Prothesoft® does not require any special milling or grinding tools, but can be machined and polished with common denture instruments.
4. The dentist is advised to place the relined denture in a water bath for approx. 30-60 minutes prior to insertion.
5. Prothesoft® softens only on contact with warm water or saliva. Otherwise, Prothesoft® remains hard like normal plastic.

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