Einleitung

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Soft relining in total prosthetics

Introduction

Based on scientific research, Prof. Dr. Nicolas Abou Tara was able to prove that Prothesoft® Forte and its base substances—as a softening relining material for dental prostheses—have neither teratogenic nor toxic properties. Prothesoft® Forte is now the third generation of this class of materials. Its composition has been significantly improved to not only reduce residual monomers, but to eliminate them completely. The material is a softening and permanently soft relining material that offers considerable advantages for dentists, dental technicians, and, above all, patients.

A permanently soft relining material should have the following properties:

1. High tissue compatibility
2. Good adhesion to the base plastic of the denture
3. Permanent elasticity and dimensional stability
4. Sufficient resilience
5. Consistent fit and abrasion resistance
6. Smooth, easily polishable surface
7. Easy processing
8. Color consistency

From a scientific point of view, the special properties of the acrylate group and its excellent bonding properties with hard acrylate plastics should be emphasized, which favors the use of such materials as relining materials. Prothesoft® Forte is based on methyl methacrylate and 2-hydroxyethyl methacrylate. Thanks to the hot polymerization process, the material is completely free of residual monomers and meets all requirements for an excellent soft relining material.

Indication and clinical relevance

The correct fit of a complete denture is a constant challenge for the dentist. The aim is to ensure that the denture remains stable without moving or lifting off the jaw ridge. The retention of a complete denture depends on several factors:

1. Anatomical jaw shape
2. Accurate impression taking
3. Precise bite registration and bite height
4. Careful dental implementation from setup to polishing
5. Fitting of the denture into the underlying jaw areas
6. Muscle-functional design of the outer surfaces
7. Creation of a suction effect

Even with optimal implementation of these parameters, the anatomical factor remains a critical issue, especially in cases of severely atrophied lower jaws. Atrophy occurs six to seven times more frequently in the lower jaw than in the upper jaw. Many affected individuals often try to improve the fit with various adhesives such as powders, creams, or films. However, these substances can cause mucosal irritation, allergic reactions, and bacterial colonization without actually achieving the desired benefit.

Limitations of surgical and implantological measures

Surgical procedures such as bone augmentation, grafts, or implants are also used to improve prosthesis retention. However, these procedures are associated with the risks of pain, complications, or even long-term failure. Neuralgia, mucosal damage, or spontaneous jaw fractures can occur, especially in severely atrophied jaws.

The adhesion and cohesion principle of Prothesoft® Forte

This problem prompted the research group led by Prof. Dr. med. Nicolas Abou Tara to develop an alternative method that does not require surgical intervention. The basis for this is the mechanism of adhesion and cohesion between the denture base (acrylate) and the saliva-covered mucosa of the jaw ridge. This creates a physical attraction between the inside of the denture and the film of saliva on the mucosa—comparable to the behavior of two moistened glass plates that can be moved sideways but are almost impossible to separate vertically.

Material properties of Prothesoft®

Prothesoft® has an oral LD₅₀ of 6.7 mg/kg in CD rats and is therefore – in accordance with Liebezeit's studies – neither toxic nor teratogenic. The material can be processed conventionally. After polymerization, it is initially hard and can be milled, ground, or polished with standard dental instruments. In the mouth, however, it behaves elastically and has a microporous structure that can absorb saliva. Within four hours after insertion, slight swelling occurs with an approx. 3% increase in weight and approx. 13% increase in volume towards the inside of the denture. This swelling causes the gaps between the denture base and the jaw ridge to be filled exactly. This creates a mechanical clamping effect in the underlying areas, which enables a stable and at the same time soft fit. The combination of elastic behavior and pumping action during chewing results in continuous absorption and release of saliva. This leads to self-cleaning of the inside of the denture and reduces the risk of pressure points or inflammation.

Summary of material advantages:

- Optimal adhesion through adhesion and cohesion effects
- Precise fit and gap filling Prevention of food retention spaces
- Valve-like closure through swelling
- Pain-free use without pressure points
- Improved suction effect between denture and jaw

Advantages of Prothesoft® Forte

a) For the dental laboratory:

- Excellent adhesion to all prosthetic materials and metals
- Long-term stability beyond the usual warranty period
- Free of residual monomers
- Hard when dry after cooling, therefore easy to work with
- No bonding agent or gloss varnish required
- Easy reworking for corrections

b) For the dentist:

- Suitable for all denture wearers, even those with allergies or anomalies
- The swelling effect improves the marginal design and denture adhesion
- Faster adjustment phase compared to rigid base materials
- High patient satisfaction – the best long-term advertising for your practice

Relining of temporary dentures has also proven itself in two-phase implantology: It prevents pressure loads and mucosal perforations that could lead to peri-implantitis.

Conclusion

Prothesoft® Forte makes the work of dentists and dental technicians considerably easier while reducing patient discomfort. The material combines proven research results with modern polymer technology developments, thus contributing significantly to the optimization of total prosthetics. A secure prosthesis fit, freedom from pain, and a long service life are the basis for successful prosthetic care—and for trust between the practitioner and the patient.